



TOPSCIENCE 众山



STERILEHYA

**Sterile
Sodium Hyaluronate**

**Master Files Registration Number:
M2025078-000**

SHANDONG TOPSCIENCE BIOTECH CO., LTD.

Address: #98 Lanshan West Road, Lanshan District,
Rizhao, Shandong, China

Tel: +86-532-80821298

E-mail: sales@topscience.cn

Website: www.topscience.cn

QINGDAO TRONGEN BIOTECH CO., LTD.

Address: #330 Songling Road, Laoshan District,
Qingdao, Shandong, China

Tel: +86-532-80821298

E-mail: sales@trongenbio.com

Website: www.topscience.cn



PRODUCT ADVANTAGES



Topsience adopts sterile production processes, eliminating the need for final sterilization steps, and ensuring product quality.



Ultra-low levels of impurities, bacterial endotoxins, residual elemental impurities, and microbial control.



Frost & Sullivan awarded Topsience Biotech Sterile Sodium Hyaluronate the honorary title of pioneer.

PRODUCTION ADVANTAGES

Completely Closed
Production Process



Tightly Controlled
Production Environment



Advanced and Complete
Production Equipments



Excellent Production
Process Standards



APPLICATIONS

Medicinal field:

Suitable for viscoelastic implantation, joint surgery, and other fillers, and can be used to prevent postoperative adhesions, etc.

Medical aesthetic field:

suitable for Class II and III medical devices and other fields.

PACKAGE AND STORAGE

1. 100g/bottle, pharmaceutical aluminum bottle with outer pharmaceutical low-density polyethylene (LDPE) bag.
2. In an airtight container, protected from light and humidity, 2-8°C.